

SPECIFICATIONS	
kV Range	40kV - 90kV
mAs Range	0.1mAs - 40mAs
Peak Power	1600W
Focal Spot	0.8mm
Input Power	Battery, Rechargeable
Collimator	In-built with high brightness lamp
Collimator Field	Max 14" X 14" @ 100cm SID
Weight	<9kG
Dimension	159 X 172 X 310 mm

Specifications are subject to changes.



OUT OF BOX CONTENTS

- X-Ray Unit Smart 1600
- 2 x Battery
- Battery Charger
- Battery Charging Dock
- 3-Meter-long Exposure Switch
- Carrying soft case

OPTIONAL ACCESSORIES

- Portable Stand
- Carrying hard case

